



Mycoplasma Detection Kits

Minerva Biolabs offers several rapid assays and solutions to detect cell culture contaminations.



>Features	Conventional PCR	qPCR	Polymerase included	Internal Control	Process Control	Positive Control	EP 2.6.7/JP G3 conformity
Venor®GeM Classic	✓	✗	✓ ✗	✓	✓	✓	✓
Venor®GeM Advance	✓	✗	✓	✓	✗	✓	✗
Venor®GeM OneStep	✓	✗	✓	✓	✗	✓	✗
Venor®GeM qEP	✗	✓	✓	✓	✓	✓	✓
Venor®GeM qOneStep	✗	✓	✓	✓	✗	✓	✗

Venor®GeM Classic



The classical version for endpoint PCR and gel electrophoresis. 20 years of experience, always improved and at the state-of-the-art of DNA-based mycoplasma detection – the optimal choice for money-savers!

Ordering information/package sizes

Without Polymerase

Cat. No. 11-1025 > 25 Reactions
Cat. No. 11-1100 > 100 Reactions
Cat. No. 11-1250 > 250 Reactions

With Polymerase

Cat. No. 11-1025G > 25 Reactions
Cat. No. 11-1100G > 100 Reactions
Cat. No. 11-1250G > 250 Reactions

Venor®GeM Advance



The pre-casted version of Venor®GeM is ready-for-gel. Just add the rehydration buffer and the sample into the pre-casted PCR tubes. After cycling the PCR products can be loaded directly onto a gel, as marker and loading buffer are already included. The best version for time-savers!

Cat. No. 11-7024 > 24 Reactions
Cat. No. 11-7048 > 48 Reactions
Cat. No. 11-7096 > 96 Reactions
Cat. No. 11-7240 > 240 Reactions

Venor®GeM OneStep



Venor®GeM OneStep is all-in-one! Do you need a robust and reliable work-horse? The Venor®GeM OneStep is the best choice as all reagents including the polymerase are freeze-dried and only need to be reactivated by adding the rehydration buffer. Best value for money!

Cat. No. 11-8025 > 25 Reactions
Cat. No. 11-8050 > 50 Reactions
Cat. No. 11-8100 > 100 Reactions
Cat. No. 11-8250 > 250 Reactions

Venor®GeM qEP



Venor®GeM qEP can be performed with any type of real-time PCR cycler able to detect the fluorescence dyes FAM™ and HEX™. Two protocols allow the fast screening of research samples or highly sensitive and robust testing of regulated biopharmaceutical samples.


Cat. No. 11-9025 > 25 Reactions
Cat. No. 11-9100 > 100 Reactions
Cat. No. 11-9250 > 250 Reactions

Venor®GeM qOneStep



TaqMan®-based qPCR assay with FAM™ and HEX™ labeled probes. Recommended for fast and reliable direct screening of cell culture supernatants. Highly sensitive and robust assay. Ready-to-use kit. Best offer!

Cat. No. 11-91025 > 25 Reactions
Cat. No. 11-91100 > 100 Reactions
Cat. No. 11-91250 > 250 Reactions



- @minervabiolabs
- Minerva Biolabs
- Minerva Biolabs
- minerva_biolabs

- Telephone
+49 (0)30 2000437-0
- E-mail
info@minerva-biolabs.com
- Web
www.minerva-biolabs.com

*FAM and HEX are trademarks of Applied Biosystems. TaqMan is a registered trademark of Roche Mol.